

Density Ball 10 ml (EN ISO 2811-2)

200011009

Product Description

Density Ball (plummet) 10 ml \pm 0.005 ml

Also known as 'Gamma Kugel'.

Method for determining the density of low viscous liquids based on the Archimedes principle. According to EN ISO 2811-2.

A body immersed in a liquid produces an upward force from which the density of the liquid can be calculated on the basis of the known volume of the body. The upward force can be read from a precise electrical balance.

Our density balls have been developed and manufactured by Solvica conform ISO 9001:2015.



Standards

EN ISO 2811-2 (formerly known as DIN 53217 part 3)

Technical Specification

Maintenance advice:

- Our density ball is precision-machined. Never drop it or knock it over
- Ensure proper cleaning this instrument after use with a suitable solvent which leaves no residues to avoid measuring errors
- Don't use any hard, abrasive materials to clean the instrument. Scratches and dents will result in permanent damage
- Please store the instrument in its original case when not in use
- Always hold the plummet in vertical position to avoid stress on the neck section which may cause damage

Density Ball 10 ml (EN ISO 2811-2)

200011009

Main Technical Parameters

Model Parameters	Density Ball (plummet) 10 ml	Density Ball (plummet) 100 ml
	Standards	EN ISO 2811-2
Material Density Ball	Premium Inox 1.4305	
Material Shaft	Anodized aluminium	
Diameter Ball	26 mm	57 mm
Volume	10 ml	100 ml
Diameter neck	1 mm	3 mm
Overall Length	260 mm	260 mm
Working Temperature	23°C	23°C

Not included (but can be ordered together with the density ball) :

- Stand for density ball (article no. 200011012)
- Electrical balance

Ordering information: article no. 2000011009 - density ball 10 ml

Density Ball 10 ml (EN ISO 2811-2)

200011009

Accessoires

- 200011012 - Stand for Density ball 10 / 100 ml

Disclaimer

The information given in this sheet is not intended to be exhaustive and any person using the product for any purpose other than that specifically recommended in this sheet without first obtaining written confirmation from us as to the suitability of the product for the intended purpose does so at his own risk. Whilst we endeavour to ensure that all advice we give about the product (whether in this sheet or otherwise) is correct we have no control over either the quality or condition of the product or the many factors affecting the use and application of the product. Therefore, unless we specifically agree in writing to do so, we do not accept any liability whatsoever or howsoever arising for the performance of the product or for any loss or damage (other than death or personal injury resulting from our negligence) arising out of the use of the product. The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development