

Muller Laboratory Grinder: 380V - 0.75 kW

BGD 650

Product Description


BGD 650 Muller Laboratory Grinder is a sturdy machine, enabling optimal and perfectly reproducible grinding and dispersion to be obtained in the laboratory. It has been designed to grind and discharge small amounts of raw materials for paint and varnish.

This unit consists of a steel frame holding two ground-glass plates, the lower of which is motor driven. A LCD screen can help operator set the total number of revolutions and automatically stops the motor once the user's defined number of revolutions has been met. Furthermore, it also can show the total after each test cycle.

An amount of raw material is spread over a limited area of the lower glass plate. The mixture is ground in successive stages, for example 50 revolutions each, at a pressure of 445N.

Between each stage, the mixture should be redistributed within the limited area.

The BGD 650 Muller Laboratory Grinder is suitable for preparing dispersions for testing mass color and tinting strength and preparing small samples for use in the quality control of pigments



Standards

ASTM D 387, ASTM D 332-B, ISO 8780-5, ISO 787/16

Technical Specification

- Controlled intelligently by microcomputer, counting accurately, reliable working status
- Touch screen and user-friendly operation interface allow operator easily to set test parameters
- Use coupling with flexible connection method to connect reducer with motor, has big torque, no step to lose and self-protection
- Grinding glass are made of heat-resisting, pressure-resisting and colorless sight glass after precise grind, with high flatness, strong and durable and no any influence to light colored sample.

Muller Laboratory Grinder: 380V - 0.75 kW

BGD 650

Main Technical Parameters

- Motor Power: 0.75KW
- Power Voltage: 380V/50Hz
- Working Speed: 75 ±5r/min
- Load: 43 kg; 64 kg; 100 kg
- Weights: 1.75 kg; 5 kg; 6.4 kg
- Grind Plate Diameter: 240 mm
- Digital Displaying Range: 0 ~ 9,999
- Overall Dimensions: 530×360×480 mm (length x width x height)
- Weight: 80 kg
- Ordering Information: BGD 650 - Muller Laboratory Grinder

Disclaimer

The information given in this sheet is not intended to be exhaustive and any person using the product for any purpose other than that specifically recommended in this sheet without first obtaining written confirmation from us as to the suitability of the product for the intended purpose does so at his own risk. Whilst we endeavour to ensure that all advice we give about the product (whether in this sheet or otherwise) is correct we have no control over either the quality or condition of the product or the many factors affecting the use and application of the product. Therefore, unless we specifically agree in writing to do so, we do not accept any liability whatsoever or howsoever arising for the performance of the product or for any loss or damage (other than death or personal injury resulting from our negligence) arising out of the use of the product. The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development