


Glass Panel Holder for BGD 2602 Glass Strips

BGD 1481

Product Description

Stainless steel glass panel holder for glass strips used as accessory for our BGD 261 and BGD 262 Drying Time Recorders



Technical Specification

This glass panel holder is used in combination with BGD 2602 Special glass strips and BGD 203 Cube Applicators.

This holding device will assist the laboratory technician to apply the paint safely and evenly on the glass strips. The holder will cover the the sharp edges of the glass strips and ensure safe application of the paint.

Glass Panel Holder for BGD 2602 Glass Strips

BGD 1481

Main Technical Parameters

Ordering information: BGD 1481 - Glass panel holder

Disclaimer

The information given in this sheet is not intended to be exhaustive and any person using the product for any purpose other than that specifically recommended in this sheet without first obtaining written confirmation from us as to the suitability of the product for the intended purpose does so at his own risk. Whilst we endeavour to ensure that all advice we give about the product (whether in this sheet or otherwise) is correct we have no control over either the quality or condition of the product or the many factors affecting the use and application of the product. Therefore, unless we specifically agree in writing to do so, we do not accept any liability whatsoever or howsoever arising for the performance of the product or for any loss or damage (other than death or personal injury resulting from our negligence) arising out of the use of the product. The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development